

## Joint Reconstruction of Phase-Cycled Balanced SSFP with Constrained Parallel Imaging

B Bilgic<sup>1,2</sup>, T Witzel<sup>1,2</sup>, H Bhat<sup>3</sup>, LL Wald<sup>1,2</sup>, K Setsompop<sup>1,2</sup>

1 Martinos Center, Charlestown, MA, USA 2 Harvard Medical School, Boston, MA, USA 3 Siemens Medical Solutions, Charlestown, MA, United States


#### **Balanced SSFP**


bSSFP has unique T2 / T1 contrast

inherent high SNR efficiency

fast imaging time: short TE & TR


#### Phase-cycled bSSFP






- But suffers from banding artifacts due to B0 inhomogeneity
- Can be mitigated by phase-cycling:
  - multiple acquisitions with different phase increment
  - this shifts location of banding artifacts

#### Phase-cycled bSSFP





#### banding-free MIP



- But suffers from banding artifacts due to B0 inhomogeneity
- Can be mitigated by phase-cycling:
  - multiple acquisitions with different phase increment
  - this shifts location of banding artifacts
  - combine cycles with Max Intensity Projection (MIP)

# max - Phase-cycled bSSFP banding-free MIP

- But suffers from banding artifacts due to B0 inhomogeneity
- Can be mitigated by phase-cycling:
  - \* multiple acquisitions increase scan time, counteracting efficiency of bSSFP
  - ❖ Parallel Imaging and Simultaneous MultiSlice (SMS) allowed up to 4-fold acceleration in phase-cycled bSSFP [1,2]


#### Joint Recon for Phase-cycled bSSFP

We propose to jointly recon phase-cycled images

- We introduce Joint L1-SPIRIT:
  - fit SPIRiT kernels [1] jointly across coils and phase-cycles
     convert banding artifacts into useful, additional spatial encoding
  - reduction in g-factor noise amplification is > 1.8-fold relative to GRAPPA
    SNR improvement is similar to 3 averages of GRAPPA recon


#### **GRAPPA** Recon

R=2 acceleration, ignoring coil and readout axes



#### Proposed: Joint L1-SPIRiT


R=2 acceleration, ignoring coil and readout axes



+ Total Variation regularization


#### Abdominal 2D bSSFP

#### four phase-cycles, acceleration R=6




#### Neuro 2D bSSFP

#### four phase-cycles, acceleration R=6



### Neuro SMS bSSFP acquisition four phase-cycles, MultiBand = 8



**SNR gain: 2 averages of Slice GRAPPA**